

ULTRA-DEEP SEQUENCING 454™ (Roche Life Sciences)

PCR and sequencing Primers

Note 1: The fusion primers used to generate amplicon libraries are each composed of four parts fused together: adaptator (in blue), key (in red), Multiplex identifiers (MIDs, in green) and template-specific sequences (in black).

Roche has designed a set of 151 MIDs. These are described in TCB No. 005-2009.

Note 2: Reverse Transcriptase sequencing by UltraDeep Sequencing 454™ is performed in 2 fragments.

•Reverse Transcriptase:

Outer primers

MJ3: 5'-AGTAGGACCTACACCTGTCA-3' (2480 to 2499)

PL1M: 5'- CCTGCTTCTGTATTTCTGCTATTAAGTCTTTTG-3' (3514 to 3546)

Fusion primers (2 fragments)

RT1:

RT1F: 5'-CGTATCGCCTCCCTCGCGCCATCAGMIDGTACCAGTAAAATTAAGCCAGGAATG-3' (2571 to 2597)

RT1R: 5'-CTATGCGCCTTGCCAGCCCGCTCAGMIDAAGCACATTGTA CTGATATCTAATCCC-3' (2970 to 2996)

RT2:

RT2F: 5'-CGTATCGCCTCCCTCGCGCCATCAGMIDTTCAGGAAGTATACTGCATTTACCATAC-3' (2919 to 2946)

RT2R: 5'-CTATGCGCCTTGCCAGCCCGCTCAGMIDATTGACAGTCCAGCTGTCTTTT-3' (3293 to 3314)

•Protease:

Outer primers

5' prot 1: 5'-TAATTTTTTAGGGAAGATCTGGCCTTCC-3' (2082 to 2109)

3' prot 1: 5'-GCAAATACTGGAGTATTGTATGGATTTTCAGG-3' (2703 to 2734)

Fusion primers

PR1F: 5'-CGTATCGCCTCCCTCGCGCCATCAGMIDTTCAGAGCAGACCAGAGCCAACAG-3' (2136 to 2158)

PR1R: 5'-CTATGCGCCTTGCCAGCCCGCTCAGMIDTACTGGTACAGTTTCAATAGGACTAAT-3' (2553 to 2580)

•Integrase:

Outer primers

INP8: 5'- TAGTAGCCAGCTGTGATAAATGTC- 3' (4336 to 4359)

INPR8: 5'- TTCCATGTTCTAATCCTCATCCTG- 3' (5082 to 5105)

Fusion primers

IN1F: 5'-CGTATCGCCTCCCTCGCGCCATCAGMIDGAAGCCATGCATGGACAAG-3' (4371 to 4389)

IN1R: 5'-CTATGCGCCTTGCCAGCCCGCTCAGMIDCTGCCATCTGTTTTCCATARTC-3' (5037 to 5058)

PCR runs on Eppendorf Mastercycler® ep Gradient S thermocycler (times have to be adapted to the thermocycler and taq polymerases used)

•Conditions for RT and prot:

- RT-PCR Step, using RT-PCR kit Titan (Roche): Outer primers

60' at 50°C

5' at 94°C

then: (30'' at 94°C; 30'' at 55°C; 1' at 68°C) x 40 cycles

- Nested PCR step, using Q5® High-Fidelity DNA polymerase (New England Biolabs®):

Fusion primers

1' at 98°C

then (10'' at 98°C; 30'' at 54°C; 45'' at 72°C) x 40 cycles

7' at 72°C

•Conditions for integrase:

- RT-PCR Step, using RT-PCR kit Titan (Roche): Outer primers

60' at 50°C

5' at 94°C

then: (30'' at 94°C; 30'' at 55°C; 1' at 68°C) x 40 cycles

- Nested PCR step, using Q5® High-Fidelity DNA polymerase (New England Biolabs®):

Fusion primers

1' at 98°C

then (10'' at 98°C; 30'' at 55°C; 45'' at 72°C) x 40 cycles

7' at 72°C

Sequencing

Library purification/quantitation, emPCR amplification, DNA library bead enrichment and sequencing run are performed according to manufacturer recommendations.